

Fact Sheet for Recipients

LightDeck Diagnostics

LightDeck® COVID-19 Total Antibody Test

Coronavirus
Disease 2019
(COVID-19)

September 4, 2020

You are being given this Fact Sheet because your sample(s) was tested for the Coronavirus Disease 2019 (COVID-19) using the LightDeck COVID-19 Total Antibody Test.

This Fact Sheet contains information to help you understand the risks and benefits of using this test for the diagnosis of COVID-19. After reading this Fact Sheet if you have questions or would like to discuss the information provided, talk to your healthcare provider.

For the most up to date information on COVID-19 please visit the CDC Coronavirus Disease 2019 (COVID-19) webpage <https://www.cdc.gov/COVID19>

What is COVID-19

COVID-19 is caused by the SARS-CoV-2 virus. The virus, which can cause mild to severe respiratory illness, has spread globally, including the United States. The current information available to characterize the spectrum of clinical illness associated with COVID-19 suggests that symptoms include cough, shortness of breath or difficulty breathing, fever, chills, muscle pain, headache, sore throat or new loss of taste or smell.

How are people tested for COVID-19?

Two kinds of tests are currently available for COVID-19: diagnostic tests and antibody tests.

- A diagnostic test tells you if you have a current infection.
- An antibody test tells you if you had a previous infection.

What is this test?

This test is an antibody test. It will help assess if you have antibodies to the virus that causes COVID-19. An antibody test may not be able to show if you have a current infection, because it can take 1-3 weeks after infection to make antibodies.

What are the known and potential risks and benefits of the test?

Potential risks include:

- Possible discomfort or other complications that can happen during sample collection.
- Possible incorrect test result (see below for more information).

Potential benefits include:

- The results, along with other information, can help your healthcare provider make informed recommendations about your care.

What does it mean if I have a positive test result?

If you have a positive test result, it is possible that you have recent or prior COVID-19 and that you have developed an

antibody response to the virus. Your healthcare provider will work with you to determine how best to care for you based on the test results along with other factors of your medical history, and your symptoms, possible exposures, and geographic location of places you have recently traveled. There is also a chance that this test can give a positive result that is wrong (a false positive result).

Even a high-performing antibody test when used in a population without many cases of COVID-19 infection may produce as many or more false results as true results because the likelihood of finding someone who has been infected is very small. Your healthcare provider will work with you to determine the likelihood of false result.

It is not known how long antibodies to SARS-CoV-2 will remain present in the body after infection. It is not known whether having antibodies to SARS-CoV-2 will protect you from getting infected again or help reduce the severity or duration of a future COVID-19 infection.

What does it mean if I have a negative test result?

A negative test result means that antibodies to the virus that causes COVID-19 were not found in your sample. However, it is possible for this test to give a negative result that is incorrect (false negative) in some people with COVID-19. A negative result may occur if you are tested early in your illness and your body hasn't had time to produce antibodies to the infection. This means that you could possibly still have COVID-19 even though the test is negative. If this is the case, your healthcare provider will consider the test result together with all other aspects of your medical history (such as symptoms, possible exposures, and geographical location of places you have recently traveled) in deciding how to care for you.

It is important that you work with your healthcare provider to help you understand the next steps you should take.

Is this test FDA-approved or cleared?

No. This test is not yet approved or cleared by the United States FDA. When there are no FDA-approved or cleared tests available, and other criteria are met, FDA can make tests available under an emergency access mechanism called an Emergency Use Authorization (EUA). This test is distributed consistent with FDA's "Policy for Coronavirus Disease-2019 Tests During the Public Health Emergency (Revised)".

What are the approved available alternatives?

There are no approved available alternative tests. FDA has issued EUAs for other antibody tests that can be found at <https://www.fda.gov/emergency-preparedness-and-response/mcm-legal-regulatory-and-policy-framework/emergency-use-authorization#2019-ncov>.

Where can I go for updates and more information? The most up-to-date information on COVID-19 is available at the CDC General webpage: <https://www.cdc.gov/COVID19> In addition, please also contact your healthcare provider with any questions/concerns.